

Digital HPLC flow meter

Designed for determine liquid flow accurately & reliably on HPLC systems.
Supplied with calibration certificate. Ideal HPLC validation tool for IQ / OQ / PQ.
Calibration is carried out gravimetrically by using water



Specification

- HPLC Liquid Flow Meter
- Range: 0.1 to 15ml/min
- Sensor calibrated for H₂O and IPA
- Resolution 0.001
- CMOS chip technology sensor
- Real time Flow reading
- Reading up to 3 place of decimal; 1.000ml/min
- Backlit LCD Display
- 12V DC operation
- Flexible to operate
- Single button operation
- Optional: Serial Printer Interface with report printout including date time stamp in the report

Marketed, Sales & Service by :

S.V. SCIENTIFIC Pvt. Ltd.

119, Sylvestra Arcade, 1st Floor, East Park Road, 11th Cross,
Malleshwaram, BENGALURU - 560003. INDIA
Phone: 080-23567255 / 23462314 / 23446548
E Mail : info@svscientific.in / sales@svscientific.in

www.svscientific.com