



Labpro Tech®

-40°C Biomedical Freezer SVF-40V328E

This products is widely used in hospitals and laboratories. Designed to store vaccines, blood plasma and many other biological materials. Installations can be found in research institutions and clinical sites in the life science, pharma, biotech, medical and electronics sectors.



Features :



Fast pull down

Better cooling performance



Coil evaporator

Better temperature uniformity



Inner-door design

4 individual inner door



Universal castor wheels

Proven heavy-duty castors with brake



Lockable doors

Standard locks with keys



High-quality steel materials

Antibacterial and corrosion resistant



Precise temp. control

Temp. control accurate to 0.1°C



5 types of alarm

Audible and visual alarms




Variable space

Height adjustable shelves




Eco-friendly

CFC-free refrigerant and energy-saving



FSF-40V328E Pull Down Test at Ambient 32°C



Marketed, Sales & Service by :

S.V. SCIENTIFIC Pvt. Ltd.

119, Sylvestra Arcade, 1st Floor, East Park Road, 11th Cross,

Malleshwaram, BENGALURU - 560003. INDIA

Phone: 080-23567255 / 23462314 / 23446548

E Mail : info@svscientific.in / sales@svscientific.in

www.svscientific.com



-40°C Biomedical Freezer SVF-40V328E

Structures :



Performance		Refrigeration System		Alarms	
Cooling method	Direct Cooling	Refrigerant	R290	High/low temperature	Y
Defrost	Manual	Temp. controller	Microprocessor	Power failure	Y
Shelf/drawer	Shelves/5	Sensor	NTC	Sensor error	Y
Temperature range	-20°C~-40°C	Display	Digital display	Controller failure	Y
Ambient temperature	10-32	External material	SPCC sprayed steel	Door ajar	Y
Noise (db)	48	Compressor QTY	1	Power failure backup	8 hours

Electrical Data		Dimensions		Other	
Voltage frequency(V/Hz)	220V 50Hz/60Hz	Package size (W*D*H)(mm)	750*760*2030	Inner door Qty	5
Power (W)	300W	Capacity (L/cu.ft)	328/11.6	Remote alarm port	Optional
Current (A)	1.6A	Exterior size (W*D*H)(mm)	690*680*1900	Access port QTY/Dia.	1/25mm
Power consumption	4kWh/24h	Interior size(W*D*H)(mm)	508*455 1393	Internal material	SPCC sprayed steel

Marketed, Sales & Service by :

S.V. SCIENTIFIC Pvt. Ltd.

119, Sylvestra Arcade, 1st Floor, East Park Road, 11th Cross,
Malleshwaram, BENGALURU - 560003. INDIA

Phone: 080-23567255 / 23462314 / 23446548

E Mail : info@svscientific.in / sales@svscientific.in

www.svscientific.com