



Labpro Tech®

Ultrasonic Cleaner / Sonicator

Characteristics :

- Microprocessor based control
- Digital controller and display
- Auto Degassing, PSP mode.
- Pharmaceutical Industries – Ensures contamination-free cleaning of lab glassware and tools.
- Hospitals – Used for pre-sterilization cleaning of surgical and medical instruments.
- R&D Laboratories – Maintains precision by cleaning sensitive lab equipment.



Main Technical Parameters :

Model (Tank Capacity in Lit)	2	3	4.5	5	7	10	15	20	30
No. of chambers	One								
Tank Size (LXWXH) in mm	240x135x65	240x135x100	300x150x100	230x230x100	300x240x100	300x240x150	500x300x100	500x300x150	500x300x200
Enclosuresize (lxwxh) In mm	265x162x200	265x162x250	325x176x240	265x265x250	320x260x320	325x270x400	525x320x400	525x320x400	525x320x450
Tank & Basket material	SS 304								
Tank Enclosure	SS 304 with matt finish								
Ultrasonic Power	50 W	100 W	100 W	100 W	150 W	200 W	300 W	400 W	500 W
Heater Power	NA	100 W	100 W	200 W	200 W	200 W	200 W	500 W	500 W
Timer	15 min	15 min	0 - 99 min						
Temp. Controller	Digital temperature controller with setting max-80°C								
Frequency	33Khz ± 5 Khz.								
Transducers	Imported make PZT type bonded at the Bottom of the tank with weld bond technique.								
Auto Degas	Pre settled for 5 Min								
Ultrasonic Generator	Advanced latest IGBT based SMPS								
Pulse Sweep Power (PSP)	For Uniform Distribution Ultrasonic Power								
Drain	½” BSP with a ball valve								
Safety Interlocking	Safety Interlocking provided with Temp for safety of the system & Overload protection by safety fuse.								
Accessories	SS 304 lid, SS 304 Basket, Drain Valve, Mains cord & User friendly operation manual.								

Marketed, Sales & Service by :

S.V. SCIENTIFIC Pvt. Ltd.

119, Sylvestra Arcade, 1st Floor, East Park Road, 11th Cross,
Malleshwaram, BENGALURU - 560003. INDIA

Phone: 080-23567255 / 23462314 / 23446548

E Mail : info@svscientific.in / sales@svscientific.in

www.svscientific.com

Email : info@svscientific.in