



Vertical Pressure Steam Sterilizer SV-HD



Characteristics :

- ✿ Fully stainless steel structure
- ✿ Hand wheel type of quick-open door structure
- ✿ Door safety lock system
- ✿ Digital display of working status, touch of key
- ✿ Self-inflating type seal
- ✿ Over temperature&over pressure auto-protection
- ✿ Safe protection of water lacking
- ✿ Automatically shut off with beep reminding after sterilization
- ✿ Auto discharge the cool air and steam discharging automatically after sterilization
- ✿ Used for sterilization of medical instruments, medical cotton products
- ✿ Fully stainless steel SUS304/AISI 304 -3mm
- ✿ 24 months warranty
- ✿ Possible to install drying system according request

Main Technical Parameters :

Model	SV-28HD	SV-35HD	SV-50HD	SV-75HD	SV-100HD
Capacity	28L (ϕ 270*500)mm	35L (ϕ 318*450)mm	50L (ϕ 340*550)mm	75L (ϕ 400*600)mm	100L (ϕ 440*650)mm
Working pressure	0.22Mpa				
Working temperature	134°C				
Max.working pressure	0.23Mpa				
Heat average	$\leq \pm 1^\circ\text{C}$				
Timer	0~99min or 0~99hours59 min				
Temp.range	105°C~134°C				
Power	2KW/ AC220V.50Hz	2.5KW/ AC220V.50Hz	3KW/ AC220V.50Hz	4.5KW/AC220V.50Hz	
Exterior size	400*370*700mm	450*450*1010mm	510*470*1130mm	560*560*1120mm	540*560*1250mm
Packing size	440*320*750mm	570*550*1150mm	590*590*1280mm	650*630*1280mm	680*630*1370mm
G.W/N.W	26/23kg	72/56kg	88kg/68kg	100kg/80kg	110kg/85kg