



Water Bath without fluid circulation 12-liter

Water bath- Digital- 12Ltr without stirrer

Outer MS, SS Inner Tank: 28x25 x 15cm (L X W X D), Temp Range: Room Temp + 5 To 99 C (Boiling water),
Outer: MS, Dismounted Lid, Hsn code: 84198990



Technical Parameters

Model	SVWTR-12-WC
Parameter	Value
Heater Capacity	Immersion Type, 1.5 Kw
Drain Provision	NO
I/P Voltage	230V AC, 50Hz
Accuracy	+/- 0.5C
Safety Thermostat	YES
Outer Dimension	33cm w*28cm L*30cm H

Marketed, Sales & Service by :

S.V. SCIENTIFIC Pvt. Ltd.

119, Sylvestra Arcade, 1st Floor, East Park Road, 11th Cross,

Malleswaram, BENGALURU - 560003. INDIA

Phone: 080-23567255 / 23462314 / 23446548

E Mail : info@svscientific.in / sales@svscientific.in

www.svscientific.com